



NRV OrthoTech Pvt Ltd
Perfection At It's Finest...



Orthopedic Implants & Instruments Company

ACL PCL Reconstruction

Mini Plates - Standard

Interlocking Nails

Small Plates - Standard

Nails / Wires / Pins

Large Plates - Standard

DHS / DCS Plates / Hip Prosthesis

Mini Locking Plates

Cannulated Screws

Small Locking Plates

Spine Implants & Instruments

Large Locking Plates

Maxillofacial Implant

External Fixators / Tools

Innovative Reliable Quality
Excellent Performance



Siwik Button (Endobutton)

TT75 / TT76

Size : 15mm to 35mm (5mm var.)

Specification

- Available in Titanium alloy (As per ISO 5832 - 2)
- Available with Sterile & Closed Loop (TT75)
- Also Comes in Non-Sterile & Without Loop (TT76)
- Width : 4.0mm
- Length : 12.0mm



Spiked Ligament Staple

TT77

Specification

- Available in Titanium alloy (As per ISO 5832 - 2)
- Available in 8mm and 11mm.
- Sharp leg points for easier penetration into cortical bone.
- Two pronged implant for superior pull out strength.



Anchors for Rotator Cuff

TT73

Types

- Suture Anchor for shoulder, Dia.5.0mm
- Suture Anchor for shoulder, Dia.3.0mm

Specification

- Available in Titanium alloy (As per ISO 5832 - 3)



Interference Screw

TT72

Dia. : 7mm, 8mm, 9mm
Size : 20mm to 35mm (5mm var.)

Specification

- Available in Titanium alloy (As per ISO 5832 - 2)



Staple Washer

TT102

Specification

- Available in Titanium alloy (As per ISO 5832 - 2) only
- Outer Diameter: 14.0mm
- Inner Diameter : 6.5mm
- No. of Spikes : 7
- Thickness : 2.0mm



ACL PCL Instrument Set



Corporate Office

NRV OrthoTech Pvt. Ltd.

C-506, Signature 2 Business Park,
Near Sanand Circle, Sarkhej,
Ahmedabad - 382210, Gujarat, INDIA.
Ph : +91-79-66173506
E-Mail : info@nrvortho.com
www.nrvortho.com