



NRV OrthoTech Pvt Ltd
Perfection At It's Finest...



Orthopedic Implants & Instruments Company

ACL PCL Reconstruction

Interlocking Nails

Nails / Wires / Pins

DHS / DCS Plates / Hip Prosthesis

Cannulated Screws

Spine Implants & Instruments

Maxillofacial Implant

Mini Plates - Standard

Small Plates - Standard

Large Plates - Standard

Mini Locking Plates

Small Locking Plates

Large Locking Plates

External Fixators / Tools

Innovative Reliable Quality
Excellent Performance



NRV Ortho

CC screw 4.0mm

SS111 / TT94

Type : Short Thread / Full Thread
Size : 10mm to 30mm (2mm var.)
35mm to 75mm (5mm var.)

Specification

- Available in Stainless Steel AISI 316 L.
- Threaded diameter : 4.0mm
- Drill bit for theraded hole : 3.2mm
- Shaft diameter : 2.7mm
- Core diameter : 2.5mm
- Head diameter : 6.0mm
- Cannulation : 1.4mm
- Guide wire : 1.2mm



NRV Ortho

CC screw 6.5mm

SS112 / TT95

Type : Thread Length 16mm / 32mm
Size : 25mm to 110mm (5mm var.)

Specification

- Available in Titanium alloy (As per ISO 5832-3) & Stainless Steel AISI 316 L.
- Threaded diameter : 6.5mm
- Drill bit for threaded hole : 4.5mm
- Shaft diameter : 4.5mm
- Core diameter : 4.3mm
- Hexagonal Socket : 3.5mm
- Head diameter : 8.0mm
- Cannulation : 2.2mm
- Guide wire : 1.8mm



NRV Ortho

CC screw 7.0mm

SS112 / TT95

Type : Thread Length 16mm / 32mm
Size : 30mm to 110mm (5mm var.)

Specification

- Available in Titanium alloy (As per ISO 5832-3) & Stainless Steel AISI 316 L.
- Threaded diameter : 7.0mm
- Drill bit for threaded hole : 4.5mm
- Shaft diameter : 4.5mm
- Core diameter : 4.3mm
- Hexagonal Socket : 3.5mm
- Head Diameter : 8.0mm
- Cannulation : 2.2mm
- Guide wire : 1.8mm



CC Washer 4.0mm & 6.5mm

SS122 / TT101

Types

- For 3.0mm Cannulated screws (6.5/3.3mm)
- For small screws (7.0/3.6mm)
- For small screws (4.5/2.2mm)
- For 4.5mm Cannulated screws (10.0/4.6mm)
- For 4.5mm to 7.3mm screws (13.0/6.6mm)
- For Locking compression plates (3.5mm)

Specification

- Available in Titanium alloy (As per ISO 5832-3) & Stainless Steel AISI 316 L.



CC Screw Box 4.0mm

1	CANULATEAD CANCELLOUS SCREW - C.C. SCREW 4mm
2	Cannulated Drill Bit
3	Cannulated Tap
4	Screw Driver for C.C.
5	Guid Wire (Simple)
6	Threaded Guide Wire
7	Double Drill Guide Sleeve
8	Direct Measuring Device
9	Box - for the Canulated Screw System



CC Screw Box 6.5mm

1	CANULATEAD CANCELLOUS SCREW - 16mm Th. 6.5mm
2	CANULATEAD CAN+CELLOUS SCREW - 32mm Th. 6.5mm
3	Cannulated Drill Bit
4	Cannulated Tap
5	Screw Driver for C.C.
6	Guide Wire (Simple)
7	Threaded Guide Wire
8	Double Drill Guide Sleeve
9	Direct Measuring Device
10	Box - for the Canulated Screw System



Herbert Screw 3.0mm

TT74

Size: 10 to 40mm (2mm var.)

Specification

Available in Titanium alloy (As per ISO 5832-2) only.



Herbert Screw Box

No.	Particular	Qty.
1	Direct Measuring Device	1
2	Double Drill Sleeve \varnothing 1.1 mm / 2.0	1
3	Drill bit \varnothing 2.0 mm, Cannulated	2
4	Guide Wire, Threaded \varnothing 1.1 mm X 150	1
5	Guide Wire, Simple \varnothing 1.1 mm X 150	1
6	Screw Driver, Cannulated	2
7	Trapezium Drill, Cannulated	1
8	Box - for Herbert Instruments	1



Corporate Office

NRV OrthoTech Pvt. Ltd.

C-506, Signature 2 Business Park,
Near Sanand Circle, Sarkhej,
Ahmedabad - 382210, Gujarat, INDIA.
Ph : +91-79-66173506
E-Mail : info@nrvortho.com
www.nrvortho.com