



**NRV OrthoTech Pvt Ltd**  
Perfection At It's Finest...



## Orthopedic Implants & Instruments Company

ACL PCL Reconstruction

Interlocking Nails

Nails / Wires / Pins

DHS / DCS Plates / Hip Prosthesis

Cannulated Screws

Spine Implants & Instruments

Maxillofacial Implant

Mini Plates - Standard

Small Plates - Standard

Large Plates - Standard

Mini Locking Plates

Small Locking Plates

Large Locking Plates

External Fixators / Tools

**Innovative Reliable Quality**  
Excellent Performance



### Dynamic Compression Plate, 2.0mm

#### Specification

- Available in Titanium alloy (As per ISO 5832-2) & Stainless Steel AISI 316 L.
- Thickness : 1.0 mm (up to 6 holes), 1.5 mm for (7 to 10 holes)
- Width : 5.0mm
- Hole spacing : 5.0mm



### Dynamic Compression Plate, 2.7mm, small

#### Specification

- Available in Stainless Steel AISI 316 L.
- Fixation with 2.7 mm screws
- Thickness : 2.0, 2.2, 2.5 & 2.7 mm
- Width : 8.0mm
- Hole spacing : 8.0mm



### T - Plate, 2.7mm

#### Specification

- Available in Stainless Steel AISI 316 L.
- Fixation with 2.7mm Screws.
- Thickness : 1.2mm
- Width : 7.0mm



### L - Plate (L&R), 2.0mm

#### Specification

- Available in Stainless Steel AISI 316 L.
- Fixation with 2.0mm screws
- Thickness : 1.0mm
- Width : 5.0mm



### Cortex Screw, 2.7mm

SS101 / TT81

Size : 6mm to 40mm (2mm var.)

#### Specification

- Available in Stainless Steel AISI 316 L.
- Thread Diameter : 2.7mm
- Drill bit for threaded role : 2.0mm
- Core diameter : 1.9mm
- Hexagonal socket : 2.5mm
- Head diameter : 5.0mm

**Corporate Office**

**NRV OrthoTech Pvt. Ltd.**

C-506, Signature 2 Business Park,  
Near Sanand Circle, Sarkhej,  
Ahmedabad - 382210, Gujarat, INDIA.  
Ph : +91-79-66173506  
E-Mail : [info@nrvortho.com](mailto:info@nrvortho.com)  
[www.nrvortho.com](http://www.nrvortho.com)